

TABLE TOP ETHYLENE OXIDE STERILIZERS MASA ÜSTÜ ETİLEN OKSİT STERİLİZATÖRÜ



AX 60 / 60 Lt.



AX 135 / 135 Lt.



**Sterilize with confidence by
Superior of Warm Sterilization at 54°C**

Cardiovascular Purpose
Kardiyovasküler



EO Dedector

Axis offers a latest technology table top warm gas diffusion sterilization system are ideal for small office-based and medium sterilization requirements. Putting what ever you wish to be sterilized in to the Liner bag along with the Sterilizing Ampoule, Humidity Chip and close with string. When Liner Bag put in the sterilizer and snap the ampoule, the sterilization is begin. The sterilized materials are useable after sterilization and aeration period. It works with only 5 and 11 gr of ampoules. Simple to Install and Operate.

Kit bazlı gaz difüzyon teknolojisi ile yeni nesil gaz siterilizatörleridir. Sterilize edilecek malzemeler, pratik olarak sterilizasyon ampul ve nem çipi ile birlikte Liner Bag içerisine konular ve ağzı kapatılır. Liner bag siterilizatörün içine yerleştirildikten sonra sterilizasyon ampülün başlığına basılarak sterilizasyon başlatılır. Sterilizasyon ve havalandırma işlemi sonrası malzemeleriniz kullanıma hazırdır. Oldukça düşük miktarlarda, 5 ve 11gr'lık EO ampülleri ile çalışır. Montaj kullanım, pratik ve kolaydır.

MODEL	VOLUME KAPASİTE	İN / İÇ			OUT / DIŞ		
		Width En	Length Boy	Depth Derinlik	Width En	Length Boy	Depth Derinlik
AX 60	60 LT	35 cm	35 cm	49 cm	37 cm	41 cm	61 cm
AX 135	135 LT	35 cm	35 cm	110 cm	37 cm	41 cm	122 cm

Each kit includes Axis sterilizing gas ampoules, Liner Bag, and Humidity chip.

Her bir kit içerisinde, Sterilizasyon ampülü, Liner Bag ve nem çipi bulunur.

Two years warranty.
İki yıl garantili.

KORDON TIP SAĞLIK ARAÇLARI LTD. ŞTİ.

Tel. 0 (232) 348 60 60 (pbx)

Fax. 0 (232) 348 60 90

Web: www.kordontip.com - www.axis.com.tr