

# AXIS STERILIZATION REEL HANGER & CUTTER

*Sterilization Packing Products*



PRODUCT CODE : AX-1300

Axis Hanger and Cutter has two separate sections. It is suitable for all kinds of Sterilization Reels. Can be used as table top or hanged on wall.

## Axis Hanger and Cutter;

- Made of Stainless Steel.
- Easy to use
- Strong durability cutting system do not generate paper dust
- Almost lifetime rotary knife
- Knife sharpens automatically while working
- Product size: 60\*38\*45 cm